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Request form for DNA typing of dogs for parentage control

Hundägare ifylles av ägaren / to be completed by the owner

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Skicka provet till SLU snarast möjligt efter provtagning. *Send the sample to SLU as soon as possible after sampling.*

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Information to the sampling veterinarian

Both the dog owner and the veterinarian should complete the form. The dog's identity must be verified at the time of sampling. Take samples of all possible parents (including mother). All puppies of the litter should be DNA typed.

Instructions for blood sampling

A blood sample (1-3 ml) is taken by a veterinarian in one EDTA-tube, clearly labelled with the dog's name, registration number or ID number. Please, turn the tube 2-3 times immediately after sampling to avoid clogging. The samples should be sent the same day or refrigerated upon sending.

Instructions for swab sampling

There are two swabs for each individual (ordered from the Animal Genetics Laboratory). Dog/puppy should not have suckled or eaten, later than 30-60 minutes before sampling. If possible, keep the puppies separated before sampling. Open the package and remove the swab making sure not to touch bristles. To collect the cells, the swab is rubbed on the inner cheek inside the mouth. Roll the swab about 10-20 times over a large area and with some pressure. Allow the swab to air dry for about 10-15 minutes at room temperature. It is important that the swab is completely dry. Place the dried swab into the provided plastic tube and close the tube with the lid. Mark the plastic tube with the dog/puppy ID (name, chip number or barcode). Repeat the procedure with the second swab on the inside of the dog/puppy's other cheek. The samples should be sent the same day or refrigerated upon sending.

Payment

Payment according to the price list is made in advance. In the message, enter the reference "670Dog" and the dog's name. Name of account holder: Sveriges Lantbruksuniversitet, Account number (IBAN): SE92 1200 0000 0128 1011 9104, BIC/Swift: DABASESX, Bank address: Norrmalmstorg 1, 111 46 Stockholm, Sweden.

Normal turnaround time is about 10-15 days after sample arrival. The results are sent by regular mail to the owner or stated recipient.

Note! Payment should be made prior to any service rendered and goods delivered.

Send the samples to the address below:

Mark the envelope: Dog parentage

Animal Genetics Laboratory
SLU
Box 7023
SE-750 07 Uppsala
Sweden