

## How to transfer data from the VR headset to your computer

This is a step-by-step guide on how you can transfer user data from the NATUREACH-app on your VR headset to your computer. To transfer all the data collected by the NATUREACH-app follow the steps below:

1. Connect your Pico 4 VR headset to your computer with a USB-C cable. You should preferably use the cable that came with the VR headset, but any USB-C cable should work.
2. Open the File Explorer on your computer.
3. Navigate to the following directory: This PC\PICO 4 Enterprise\Internal shared storage\Android\data\com.VirtualDawn.Natureach360\files
4. There you will see a folder called Data. Open the folder.
5. The folder contains .csv files containing the user data. The file names indicate the date and time the data was collected. Select the files you want to transfer.
6. Once you have selected all the files you want to transfer right-click one of the files. Select Copy.
7. Navigate to the preferred destination on your computer, where you want to copy the data to.
8. Transfer the data into the folder by right-clicking and selecting Paste.